

## Product datasheet for **TA804828M**

### Calpain 5 (CAPN5) Mouse Monoclonal Antibody [Clone ID: OTI1F10]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI1F10
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 274-578 of human CAPN5 (NP_004046) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	73 kDa
Gene Name:	calpain 5
Database Link:	<a href="#">NP_004046</a> <a href="#">Entrez Gene 12337 Mouse</a> <a href="#">Entrez Gene 171495 Rat</a> <a href="#">Entrez Gene 726 Human</a> <a href="#">O15484</a>
Background:	Calpains are calcium-dependent cysteine proteases involved in signal transduction in a variety of cellular processes. A functional calpain protein consists of an invariant small subunit and 1 of a family of large subunits. CAPN5 is one of the large subunits. Unlike some of the calpains, CAPN5 and CAPN6 lack a calmodulin-like domain IV. Because of the significant similarity to <i>Caenorhabditis elegans</i> sex determination gene tra-3, CAPN5 is also called as HTRA3. [provided by RefSeq, Jul 2008]


[View online »](#)

**Synonyms:** ADNIV; HTRA3; nCL-3; VRNI

**Protein Families:** Druggable Genome, Protease

**Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY CAPN5 (Cat# [RC202045], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CAPN5 (Cat# [TA804828]). Positive lysates [LY418247] (100ug) and [LC418247] (20ug) can be purchased separately from OriGene.