

## Product datasheet for **TA805137BM**


### FLIP (CFLAR) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1B1]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1B1
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 1-272 of human CFLAR (NP_003870) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	55.2 kDa
Gene Name:	CASP8 and FADD like apoptosis regulator
Database Link:	<a href="#">NP_003870</a> <a href="#">Entrez Gene 8837 Human</a> <a href="#">O15519</a>
Synonyms:	c-FLIP; c-FLIPL; c-FLIPR; c-FLIPS; CASH; CASP8AP1; Casper; CLARP; FLAME; FLAME-1; FLAME1; FLIP
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Apoptosis



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CFLAR ([RC200608], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CFLAR. Positive lysates [LY418376] (100ug) and [LC418376] (20ug) can be purchased separately from OriGene.