

## Product datasheet for **TA805300BM**


### **TWIST2 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI6G9]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI6G9
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human TWIST2 (NP_476527) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	twist family bHLH transcription factor 2
Database Link:	<a href="#">NP_476527</a> <a href="#">Entrez Gene 13345 Mouse</a> <a href="#">Entrez Gene 59327 Rat</a> <a href="#">Entrez Gene 117581 Human</a> <a href="#">Q8WVJ9</a>
Synonyms:	AMS; BBRSAY; bHLHa39; DERMO1; FFDD3; SETLSS
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TWIST2 ([RC205006], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TWIST2. Positive lysates [LY403298] (100ug) and [LC403298] (20ug) can be purchased separately from OriGene.