

Product datasheet for **TA805346M**

HIP55 (DBNL) Mouse Monoclonal Antibody [Clone ID: OTI7B7]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI7B7
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human DBNL (NP_054782) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	48.1 kDa
Gene Name:	drebrin like
Database Link:	NP_054782 Entrez Gene 13169 Mouse Entrez Gene 83527 Rat Entrez Gene 28988 Human Q9UJU6
Synonyms:	ABP1; HIP-55; HIP55; SH3P7



[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY DBNL (Cat# [RC203435], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DBNL (Cat# [TA805346]). Positive lysates [LY402274] (100ug) and [LC402274] (20ug) can be purchased separately from OriGene.