

Product datasheet for TA805928

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Factor I (CFI) Mouse Monoclonal Antibody [Clone ID: OTI7C8]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI7C8

Applications: IHC

Recommended Dilution: IHC 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 340-583 of human

CFI(NP_000195) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 63.4 kDa

Gene Name: complement factor I

Database Link: NP 000195

Entrez Gene 3426 Human

P05156



Background:


This gene encodes a serine proteinase that is essential for regulating the complement cascade. The encoded preproprotein is cleaved to produce both heavy and light chains, which are linked by disulfide bonds to form a heterodimeric glycoprotein. This heterodimer can cleave and inactivate the complement components C4b and C3b, and it prevents the assembly of the C3 and C5 convertase enzymes. Defects in this gene cause complement factor I deficiency, an autosomal recessive disease associated with a susceptibility to pyogenic infections. Mutations in this gene have been associated with a predisposition to atypical hemolytic uraemic syndrome, a disease characterized by acute renal failure, microangiopathic hemolytic anemia and thrombocytopenia. Primary glomerulonephritis with immmune deposits is another condition associated with mutation of this gene. [provided by RefSeq, Jul 2008]

Synonyms: AHUS3; ARMD13; C3b-INA; C3BINA; FI; IF; KAF

Protein Families: Druggable Genome, Protease, Secreted Protein


Protein Pathways: Complement and coagulation cascades

Product images:




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-CFI mouse monoclonal antibody. This figure shows negative staining. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA805928) (1:2000)





Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-CFI mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA805928) (1:2000)



Immunohistochemical staining of paraffinembedded Human endometrium tissue within the normal limits using anti-CFI mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA805928) (1:2000)