

## Product datasheet for **TA806069AM**


### GRK4 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI9A3]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI9A3
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 1-245 of human GRK4(NP_892027) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	66.4 kDa
Gene Name:	G protein-coupled receptor kinase 4
Database Link:	<a href="#">NP_892027</a> <a href="#">Entrez Gene 2868 Human</a> <a href="#">P32298</a>
Synonyms:	GPRK2L; GPRK4; GRK4a; IT11
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Chemokine signaling pathway, Endocytosis



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY GRK4 (Cat# [RC216602], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GRK4(Cat# [TA806069]). Positive lysates [LY405296] (100ug) and [LC405296] (20ug) can be purchased separately from OriGene.