

## Product datasheet for **TA806154S**


### Steroidogenic Factor 1 (NR5A1) Mouse Monoclonal Antibody [Clone ID: OTI10E2]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI10E2
Applications:	IHC
Recommended Dilution:	IHC 1:200
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 1-262 of human NR5A1(NP_004950) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	51.5 kDa
Gene Name:	nuclear receptor subfamily 5 group A member 1
Database Link:	<a href="#">NP_004950</a> <a href="#">Entrez Gene 26423 Mouse</a> <a href="#">Entrez Gene 83826 Rat</a> <a href="#">Entrez Gene 2516 Human</a> <a href="#">Q13285</a>
Synonyms:	AD4BP; ELP; FTZ1; FTZF1; hSF-1; POF7; SF-1; SF1; SPGF8; SRXY3
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**Product images:**

Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human endometrium tissue from three different patients using anti-NR5A1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3min, [TA806154]) (1:200)



Immunohistochemical staining of paraffin-embedded cervical carcinoma tissue using anti-NR5A1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3min, [TA806154]) (1:200)