

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TA806249BM

## PCBP1 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI3B9]

## **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI3B9
Applications:	IHC, WB
Recommended Dilution:	WB 1:2000, IHC 1:150
Reactivity:	Human, Mouse, Rat
Host:	Mouse
lsotype:	lgG2b
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 1-225 of human PCBP1(NP_006187) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	37.3 kDa
Gene Name:	poly(rC) binding protein 1
Database Link:	<u>NP_006187</u> <u>Entrez Gene 23983 MouseEntrez Gene 500242 RatEntrez Gene 5093 Human</u> <u>Q15365</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

### CRIGENE PCBP1 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI3B9] – TA806249BM


Background:This intronless gene is thought to have been generated by retrotransposition of a fully<br/>processed PCBP-2 mRNA. This gene and PCBP-2 have paralogues (PCBP3 and PCBP4) which<br/>are thought to have arisen as a result of duplication events of entire genes. The protein<br/>encoded by this gene appears to be multifunctional. It along with PCBP-2 and hnRNPK<br/>corresponds to the major cellular poly(rC)-binding protein. It contains three K-homologous<br/>(KH) domains which may be involved in RNA binding. This encoded protein together with<br/>PCBP-2 also functions as translational coactivators of poliovirus RNA via a sequence-specific<br/>interaction with stem-loop IV of the IRES and promote poliovirus RNA replication by binding<br/>to its 5'-terminal cloverleaf structure. It has also been implicated in translational control of the<br/>15-lipoxygenase mRNA, human Papillomavirus type 16 L2 mRNA, and hepatitis A virus RNA.<br/>The encoded protein is also suggested to play a part in formation of a sequence-specific<br/>alpha-globin mRNP complex which is associated with alpha-globin mRNA stability. [provided<br/>by RefSeq, Jul 2008]

Synonyms: HEL-S-85; hnRNP-E1; hnRNP-X; HNRPE1; HNRPX

Spliceosome


Protein Pathways:

### **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PCBP1 ([RC207878], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PCBP1. Positive lysates [LY401867] (100ug) and [LC401867] (20ug) can be purchased separately from OriGene.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Immunohistochemical staining of paraffinembedded Human colon tissue within the normal limits using anti-PCBP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, [TA806249])

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US