

Product datasheet for TA806588

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

NR1D2 Mouse Monoclonal Antibody [Clone ID: OTI3D4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3D4

Applications: WB

Recommended Dilution: WB 1:2000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 1-294 of human

NR1D2(NP_005117) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: nuclear receptor subfamily 1 group D member 2

Database Link: NP 005117

Entrez Gene 259241 RatEntrez Gene 353187 MouseEntrez Gene 9975 Human


Q14995

Synonyms: BD73; EAR-1R; RVR

Protein Families: Druggable Genome, Nuclear Hormone Receptor, Transcription Factors



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY NR1D2 (Cat# [RC206338], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NR1D2(Cat# TA806588). Positive lysates [LY417483] (100ug) and [LC417483] (20ug) can be purchased separately from OriGene.