

Product datasheet for **TA806758**

CD33 Mouse Monoclonal Antibody [Clone ID: OTI7G8]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI7G8
Applications:	FC, WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human CD33 (NP_001763) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	CD33 molecule
Database Link:	NP_001763 Entrez Gene 945 Human P20138
Synonyms:	p67; SIGLEC-3; SIGLEC3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Hematopoietic cell lineage




[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CD33 ([RC207023], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CD33. Positive lysates [LY400667] (100ug) and [LC400667] (20ug) can be purchased separately from OriGene.



Flow cytometric analysis of living 293T cells transfected with CD33 overexpression plasmid ([RC207023]), Red/empty vector ([PS100001], Blue) using anti-CD33 antibody ([TA506758]). Cells incubated with a non-specific antibody (Green) were used as isotype control (1:100).