

Product datasheet for TA807198AM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CD3E Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: UCHT1]

Product data:

Product Type: Primary Antibodies

Clone Name: UCHT1
Applications: FC, WB

Recommended Dilution: WB 1:2000, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human infant thymocytes and lymphocytes from a patient with Sézary disease

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Biotin

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 20.7 kDa

Gene Name: CD3e molecule

Database Link: NP 000724

Entrez Gene 916 Human

P07766





Background:


The protein encoded by this gene is the CD3-epsilon polypeptide, which together with CD3-gamma, -delta and -zeta, and the T-cell receptor alpha/beta and gamma/delta heterodimers, forms the T-cell receptor-CD3 complex. This complex plays an important role in coupling antigen recognition to several intracellular signal-transduction pathways. The genes encoding the epsilon, gamma and delta polypeptides are located in the same cluster on chromosome 11. The epsilon polypeptide plays an essential role in T-cell development. Defects in this gene cause immunodeficiency. This gene has also been linked to a susceptibility to type I diabetes in women. [provided by RefSeq, Jul 2008]

Synonyms: IMD18; T3E; TCRE


Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Hematopoietic cell lineage, Primary immunodeficiency, T cell receptor signaling pathway

Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CD3 ([RC208276], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CD3 (1:2000). Positive lysates [LY400242] (100ug) and [LC400242] (20ug) can be purchased separately from OriGene.




Flow cytometric Analysis of live Jurkat cells, using anti-CD3 antibody ([TA807198]), (Red), compared to a nonspecific negative control antibody (Mouse-IGg), (Blue) (1:100).





Flow cytometric Analysis of living RBC-lysed human peripheral blood cells, using anti-CD3 antibody ([TA807198], Red), compared to an IgG isotype control (green), and negative control (PBS, Blue) (1:100).



Flow cytometric Analysis of living peripheral blood mononuclear cells (PBMC), using anti-CD3 antibody ([TA807198], Red), compared to an IgG isotype control (Blank) (1:100).