

## Product datasheet for **TA807377AM**


### **C19orf80 (ANGPTL8) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI4F1]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI4F1
Applications:	IHC
Recommended Dilution:	IHC 1:150
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 22-198 of human C19ORF80(NP_061157) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	21.9 kDa
Gene Name:	angiopoietin like 8
Database Link:	<a href="#">NP_061157</a> <a href="#">Entrez Gene 55908 Human</a> <a href="#">Q6UXH0</a>
Synonyms:	C19orf80; PRO1185; PVPA599; RIFL; TD26



[View online »](#)

**Product images:**

Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human endometrium tissue using anti-C19ORF80 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA807377]) (1:150)