

### **Product datasheet for TA807534**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## MelanA (MLANA) Mouse Monoclonal Antibody [Clone ID: OTI2E7]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI2E7

**Applications:** IF, IHC, WB

Recommended Dilution: WB 1:2000, IHC 1:1200

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 48-118 of human

MLANA(NP\_005502) produced in E.coli.

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

Predicted Protein Size: 13 kDa

Gene Name: melan-A

Database Link: NP 005502

Entrez Gene 2315 Human

Q16655


Synonyms: MART-1; MART1

Protein Families: Transmembrane






# **Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MLANA ([RC204190], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MLANA (1:2000). Positive lysates [LY417261] (100ug) and [LC417261] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 3 different cell lines by using anti-MLANA monoclonal antibody (1:200).






Immunohistochemical staining of paraffinembedded Human melanoma tissue using anti-MLANA mouse monoclonal antibody. (Heatinduced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA807534) (1:1200)








Immunocytochemistry staining of SK-MEL-5 cells using anti-MLANA mouse monoclonal antibody (TA807534) (Left). The right is negative control (1:10000).

SK-MEL-28 MCF7





Immunocytochemistry staining of SK-MEL-28 cells using anti-MLANA mouse monoclonal antibody (TA807534) (Left). The right is negative control (1:10000).