

Product datasheet for **TA807769**

CGK2 (PRKG2) Mouse Monoclonal Antibody [Clone ID: OTI15E6]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI15E6
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 40-320 of human PRKG2 (NP_006250) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	87.3 kDa
Gene Name:	protein kinase, cGMP-dependent, type II
Database Link:	NP_006250 Entrez Gene 19092 Mouse Entrez Gene 25523 Rat Entrez Gene 5593 Human Q13237
Background:	This gene encodes a protein that belongs to the serine/threonine protein kinase family of proteins. The encoded protein plays a role in the regulation of fluid balance in the intestine. A similar protein in mouse is thought to regulate differentiation and proliferation of cells in the colon. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]
Synonyms:	cGK2; cGKII; PRKGR2




[View online »](#)

Protein Families: Druggable Genome, Protein Kinase

Protein Pathways: Gap junction, Long-term depression, Olfactory transduction

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PRKG2 ([RC215879], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PRKG2 (1:2000). Positive lysates [LY401884] (100ug) and [LC401884] (20ug) can be purchased separately from OriGene.