

## Product datasheet for **TA807783S**


### UBL3 Mouse Monoclonal Antibody [Clone ID: OTI1G4]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1G4
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 1-114 of human UBL3(NP_009037) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	13 kDa
Gene Name:	ubiquitin like 3
Database Link:	<a href="#">NP_009037</a> <a href="#">Entrez Gene 24109 Mouse</a> <a href="#">Entrez Gene 363869 Rat</a> <a href="#">Entrez Gene 5412 Human</a> <a href="#">O95164</a>
Synonyms:	HCG-1; PNSC1
Protein Families:	Druggable Genome



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY UBL3 (Cat# [RC208453], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UBL3 (Cat# [TA807783])(1:2000). Positive lysates [LY416200] (100ug) and [LC416200] (20ug) can be purchased separately from OriGene.