

Product datasheet for **TA808378AM**


ANKRA2 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI8A9]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI8A9
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ANKRA2 (NP_075526) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	34.1 kDa
Gene Name:	ankyrin repeat family A member 2
Database Link:	NP_075526 Entrez Gene 68558 Mouse Entrez Gene 294679 Rat Entrez Gene 57763 Human Q9H9E1
Synonyms:	ANKRA



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ANKRA2 (Cat# [RC204384], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ANKRA2 (Cat# [TA808378])(1:2000). Positive lysates [LY411499] (100ug) and [LC411499] (20ug) can be purchased separately from OriGene.