

## Product datasheet for **TA808408**

### **FAM104B Mouse Monoclonal Antibody [Clone ID: OTI1D10]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI1D10
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FAM104B(NP_612371) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	12.9 kDa
Gene Name:	family with sequence similarity 104 member B
Database Link:	<a href="#">NP_612371</a> <a href="#">Entrez Gene 90736 Human</a> <a href="#">Q5XKR9</a>
Synonyms:	CXorf44



[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY FAM104B (Cat# [RC205045], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FAM104B (Cat# TA808408)(1:2000). Positive lysates [LY408633] (100ug) and [LC408633] (20ug) can be purchased separately from OriGene.