


## Product datasheet for **TA808629S**

### **GAS8 Mouse Monoclonal Antibody [Clone ID: OTI10G2]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI10G2
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human GAS8 (NP_001472) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	56.2 kDa
Gene Name:	growth arrest specific 8
Database Link:	<a href="#">NP_001472</a> <a href="#">Entrez Gene 104346 Mouse</a> <a href="#">Entrez Gene 361438 Rat</a> <a href="#">Entrez Gene 2622 Human</a> <a href="#">O95995</a>
Background:	This gene includes 11 exons spanning 25 kb and maps to a region of chromosome 16 that is sometimes deleted in breast and prostate cancer. The second intron contains an apparently intronless gene, C16orf3, that is transcribed in the opposite orientation. This gene is a putative tumor suppressor gene. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2013]
Synonyms:	CILD33; DRC4; GAS11


[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY GAS8 (Cat# [RC211150], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GAS8 (Cat# [TA808629])(1:2000). Positive lysates [LY419922] (100ug) and [LC419922] (20ug) can be purchased separately from OriGene.