

Product datasheet for TA808685S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

EGLN2 Mouse Monoclonal Antibody [Clone ID: OTI5E11]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5E11

Applications: WB

Recommended Dilution: WB 1:2000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human EGLN2(NP_542770) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: egl-9 family hypoxia inducible factor 2

Database Link: NP 542770

Entrez Gene 112406 MouseEntrez Gene 308457 RatEntrez Gene 112398 Human

Q96KS0

Background: The hypoxia inducible factor (HIF) is a transcriptional complex that is involved in oxygen

homeostasis. At normal oxygen levels, the alpha subunit of HIF is targeted for degration by prolyl hydroxylation. This gene encodes an enzyme responsible for this post-translational modification. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the upstream RAB4B (RAB4B, member RAS

oncogene family) gene. [provided by RefSeq, Feb 2011]

Synonyms: EIT6; HIF-PH1; HIFPH1; HPH-1; HPH-3; PHD1






Protein Families: Druggable Genome

Protein Pathways: Pathways in cancer, Renal cell carcinoma

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY EGLN2 ([RC219191], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-EGLN2 (1:2000). Positive lysates [LY409048] (100ug) and [LC409048] (20ug) can be purchased separately from OriGene.