

Product datasheet for TA808702

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Her2 (ERBB2) Mouse Monoclonal Antibody [Clone ID: 5G12-21]

Product data:

Product Type: Primary Antibodies

Clone Name: 5G12-21

Applications: FC

Recommended Dilution: FC: 1:50

Reactivity: Human, Mouse, Rat

Host: Mouse

Isotype: IgG2b, kappa
Clonality: Monoclonal

Immunogen: HER2 extracellular domain (aa 23-652) expressed in HEK293 cells with a C-terminal DDK tag

Formulation: PBS with 50% glycerol, 1% BSA and 0.02% sodium azide

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: erb-b2 receptor tyrosine kinase 2

Database Link: NP 004439

Entrez Gene 13866 MouseEntrez Gene 24337 RatEntrez Gene 2064 Human

P04626

Synonyms: CD340; HER-2; HER-2/neu; HER2; MLN 19; NEU; NGL; TKR1

Protein Families: Druggable Genome, Protein Kinase, Transmembrane

Protein Pathways: Adherens junction, Bladder cancer, Calcium signaling pathway, Endometrial cancer, ErbB


signaling pathway, Focal adhesion, Non-small cell lung cancer, Pancreatic cancer, Pathways in

cancer. Prostate cancer





Product images:



Flow cytometric analysis of live MCF7 cells (non-fixed, non-permeabilized), using anti-HER2/ERBB2 antibody (TA808702 at 1:50, Brown), compared to IgG2b isotype control (Black).