

## Product datasheet for **TA808851S**


### Upstream Binding Protein 1 (UBP1) Mouse Monoclonal Antibody [Clone ID: OTI9G7]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI9G7
Applications:	IHC
Recommended Dilution:	IHC 1:150
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 47-269 of human UBP1(NP_055332) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	60.3 kDa
Gene Name:	upstream binding protein 1 (LBP-1a)
Database Link:	<a href="#">NP_055332</a> <a href="#">Entrez Gene 22221 Mouse</a> <a href="#">Entrez Gene 301038 Rat</a> <a href="#">Entrez Gene 7342 Human</a> <a href="#">Q9NZI7</a>
Synonyms:	LBP-1a; LBP-1B; LBP1A; LBP1B
Protein Families:	Transcription Factors




[View online »](#)

**Product images:**

Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human breast tissue using anti-UBP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA808851]) (1:150)



Immunohistochemical staining of paraffin-embedded Human thyroid tissue within the normal limits using anti-UBP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA808851]) (1:150)



Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human endometrium tissue using anti-UBP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA808851]) (1:150)