

Product datasheet for **TA808853AM**

Upstream Binding Protein 1 (UBP1) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI8F1]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI8F1
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 47-269 of human UBP1(NP_055332) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	60.3 kDa
Gene Name:	upstream binding protein 1 (LBP-1a)
Database Link:	NP_055332 Entrez Gene 22221 MouseEntrez Gene 301038 RatEntrez Gene 7342 Human Q9NZI7
Synonyms:	LBP-1a; LBP-1B; LBP1A; LBP1B
Protein Families:	Transcription Factors



[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY UBP1 (Cat# [RC208469], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UBP1 (Cat# [TA808853])(1:2000). Positive lysates [LY415219] (100ug) and [LC415219] (20ug) can be purchased separately from OriGene.