

Product datasheet for TA808947S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

DECR1 Mouse Monoclonal Antibody [Clone ID: OTI5D6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5D6

Applications: FC, IHC, WB

Recommended Dilution: WB 1:2000, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human DECR1 (NP_001350) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 32.1 kDa

Gene Name: 2,4-dienoyl-CoA reductase 1, mitochondrial

Database Link: NP 001350

Entrez Gene 1666 Human

Q16698

Background: This gene encodes an accessory enzyme which participates in the beta-oxidation and

metabolism of unsaturated fatty enoyl-CoA esters. [provided by RefSeq, Jul 2008]


Synonyms: DECR; NADPH; SDR18C1

Protein Families: Druggable Genome






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DECR1 ([RC211257], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DECR1 (1:2000). Positive lysates [LY419973] (100ug) and [LC419973] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-DECR1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA808947]) (1:150)





Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-DECR1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA808947]) (1:150)



Flow cytometric Analysis of permeabilized Hela cells, using anti-DECR1 antibody ([TA808947]), (Red), compared to negative control (PBS), (Blue) (1:100).