

### Product datasheet for TA808966BM

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **DECR1 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI16C2]**

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI16C2

Applications: WB

Recommended Dilution: WB 1:2000

Reactivity: Human

**Host:** Mouse

**Isotype:** IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human DECR1 (NP\_001350) produced in HEK293T

cell

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 32.1 kDa

**Gene Name:** 2,4-dienoyl-CoA reductase 1, mitochondrial

Database Link: NP 001350

Entrez Gene 1666 Human

016698

**Background:** This gene encodes an accessory enzyme which participates in the beta-oxidation and

metabolism of unsaturated fatty enoyl-CoA esters. [provided by RefSeq, Jul 2008]


Synonyms: DECR; NADPH; SDR18C1

**Protein Families:** Druggable Genome





# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DECR1 ([RC211257], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DECR1 (1:2000). Positive lysates [LY419973] (100ug) and [LC419973] (20ug) can be purchased separately from OriGene.