

## Product datasheet for **TA809000BM**


### **RAB3IP Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI5G4]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI5G4
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human RAB3IP (NP_071901) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	RAB3A interacting protein
Database Link:	<a href="#">NP_071901</a> <a href="#">Entrez Gene 29885 Rat</a> <a href="#">Entrez Gene 216363 Mouse</a> <a href="#">Entrez Gene 117177 Human</a> <a href="#">Q96QF0</a>
Synonyms:	RABIN3; RABIN8



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RAB3IP [RC218153], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RAB3IP (1:2000). Positive lysates [LY411671] (100ug) and [LC411671] (20ug) can be purchased separately from OriGene.