

## Product datasheet for **TA809102BM**


### ZNF230 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI7C7]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI7C7
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ZNF230 (NP_006291) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	54.4 kDa
Gene Name:	zinc finger protein 230
Database Link:	<a href="#">NP_006291</a> <a href="#">Entrez Gene 7773 Human</a> <a href="#">Q9UIEQ</a>
Background:	May be involved in transcriptional regulation. [UniProtKB/Swiss-Prot Function]
Synonyms:	FDZF2
Protein Families:	Transcription Factors



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ZNF230 (Cat# [RC208201], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF230 (Cat# [TA809102])(1:500). Positive lysates [LY416742] (100ug) and [LC416742] (20ug) can be purchased separately from OriGene.