

## Product datasheet for **TA809206S**


### UEVLD Mouse Monoclonal Antibody [Clone ID: OTI8G3]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI8G3
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 131-379 of human UEVLD(NP_060784) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	42.1 kDa
Gene Name:	UEV and lactate/malate dehydrogenase domains
Database Link:	<a href="#">NP_060784</a> <a href="#">Entrez Gene 55293 Human</a> <a href="#">Q8IX04</a>
Synonyms:	ATTP; UEV3

[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY UEVLD ([RC215907], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UEVLD (1:2000). Positive lysates [LY413157] (100ug) and [LC413157] (20ug) can be purchased separately from OriGene.