

## Product datasheet for **TA809339AM**

### **C1orf83 (TCEANC2) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI2B11]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI2B11
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human TCEANC2(NP_694580) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	24 kDa
Gene Name:	transcription elongation factor A N-terminal and central domain containing 2
Database Link:	<a href="#">NP_694580</a> <a href="#">Entrez Gene 66526 Mouse</a> <a href="#">Entrez Gene 366431 Rat</a> <a href="#">Entrez Gene 127428 Human</a> <a href="#">Q96MN5</a>
Synonyms:	C1orf83
Protein Families:	Transcription Factors



[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TCEANC2 ([RC206749], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TCEANC2 (1:2000). Positive lysates [LY407172] (100ug) and [LC407172] (20ug) can be purchased separately from OriGene.