

## Product datasheet for **TA809340BM**

### **C1orf83 (TCEANC2) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI3E3]**


#### **Product data:**

|                         |  |
|-------------------------|--|
| Product Type:           | Primary Antibodies   |
| Clone Name:             | OTI3E3   |
| Applications:           | WB   |
| Recommended Dilution:   | WB 1:2000  |
| Reactivity:             | Human, Mouse, Rat  |
| Host:                   | Mouse  |
| Isotype:                | IgG1   |
| Clonality:              | Monoclonal   |
| Immunogen:              | Full length human recombinant protein of human TCEANC2(NP_694580) produced in E.coli.  |
| Formulation:            | PBS (pH 7.3) containing 1% BSA, 50% glycerol.  |
| Concentration:          | 0.5 mg/ml  |
| Purification:           | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)  |
| Conjugation:            | HRP  |
| Storage:                | Store at -20°C as received.  |
| Stability:              | Stable for 12 months from date of receipt.   |
| Predicted Protein Size: | 24 kDa   |
| Gene Name:              | transcription elongation factor A N-terminal and central domain containing 2   |
| Database Link:          | <a href="#">NP_694580</a><br><a href="#">Entrez Gene 66526 Mouse</a> <a href="#">Entrez Gene 366431 Rat</a> <a href="#">Entrez Gene 127428 Human</a><br><a href="#">Q96MN5</a> |
| Synonyms:               | C1orf83  |
| Protein Families:       | Transcription Factors  |



[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TCEANC2 ([RC206749], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TCEANC2 (1:2000). Positive lysates [LY407172] (100ug) and [LC407172] (20ug) can be purchased separately from OriGene.