

## Product datasheet for **TA809437**


### RhoGAP (ARHGAP5) Mouse Monoclonal Antibody [Clone ID: OTI9H2]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI9H2
Applications:	WB
Recommended Dilution:	WB 1:500~2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 948-1261 of human ARHGAP5 (NP_001164) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	Rho GTPase activating protein 5
Database Link:	<a href="#">NP_001164</a> <a href="#">Entrez Gene 11855 Mouse</a> <a href="#">Entrez Gene 299012 Rat</a> <a href="#">Entrez Gene 394 Human</a> <a href="#">Q13017</a>
Background:	Rho GTPase activating protein 5 negatively regulates RHO GTPases, a family which may mediate cytoskeleton changes by stimulating the hydrolysis of bound GTP. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Synonyms:	GFI2; p190-B; p190BRhoGAP; RhoGAP5
Protein Pathways:	Focal adhesion, Leukocyte transendothelial migration



[View online »](#)

**Product images:**


HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ARHGAP5 (Cat# [RC217203], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ARHGAP5 (Cat# TA809437)(1:2000).



Western blot analysis of extracts (35ug) from 8 different cell lines by using anti-ARHGAP5 monoclonal antibody (1:500).