

## Product datasheet for **TA809517S**

### Cezanne (OTUD7B) Mouse Monoclonal Antibody [Clone ID: OTI3D10]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI3D10
Applications:	WB
Recommended Dilution:	WB 1:500~2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 449-729 of human OTUD7B (NP_064590) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	92.3 kDa
Gene Name:	OTU deubiquitinase 7B
Database Link:	<a href="#">NP_064590</a> <a href="#">Entrez Gene 229603 Mouse</a> <a href="#">Entrez Gene 310677 Rat</a> <a href="#">Entrez Gene 56957 Human</a> <a href="#">Q6GQQ9</a>
Synonyms:	CEZANNE; ZA20D1




[View online »](#)

**Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY OTUD7B [RC209932], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-OTUD7B (1:2000). Positive lysates [LY412591] (100ug) and [LC412591] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 3 different cell lines by using anti-OTUD7B monoclonal antibody (1:500).