

## Product datasheet for **TA809679BM**


### NUDT4 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1C2]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1C2
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human NUDT4 (NP_061967) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	20.1 kDa
Gene Name:	nudix hydrolase 4
Database Link:	<a href="#">NP_061967</a> <a href="#">Entrez Gene 11163 Human</a> <a href="#">A0A024RBG1</a>
Synonyms:	DIPP2; DIPP2alpha; DIPP2beta; HDCMB47P
Protein Families:	Druggable Genome



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY NUDT4 (Cat# [RC204100], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NUDT4 (Cat# [TA809679])(1:2000). Positive lysates [LY412753] (100ug) and [LC412753] (20ug) can be purchased separately from OriGene.