

Product datasheet for **TA809685**


ZMAT4 Mouse Monoclonal Antibody [Clone ID: OTI9E11]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI9E11
Applications:	WB
Recommended Dilution:	WB 1:500~2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ZMAT4 (NP_078921) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	25.7 kDa
Gene Name:	zinc finger matrin-type 4
Database Link:	NP_078921 Entrez Gene 79698 Human Q9H898
Synonyms:	FLJ13842



[View online »](#)

Product images:


HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ZMAT4 (Cat# [RC218769], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZMAT4 (Cat# TA809685)(1:2000). Positive lysates [LY403013] (100ug) and [LC403013] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from HELA cell line by using anti-ZMAT4 monoclonal antibody (1:500).