

## Product datasheet for **TA809685BM**

### ZMAT4 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI9E11]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI9E11
Applications:	WB
Recommended Dilution:	WB 1:500~2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ZMAT4 (NP_078921) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	25.7 kDa
Gene Name:	zinc finger matrin-type 4
Database Link:	<a href="#">NP_078921</a> <a href="#">Entrez Gene 79698 Human</a> <a href="#">Q9H898</a>
Synonyms:	FLJ13842




[View online »](#)

**Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ZMAT4 (Cat# [RC218769], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZMAT4 (Cat# [TA809685])(1:2000). Positive lysates [LY403013] (100ug) and [LC403013] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from HELA cell line by using anti-ZMAT4 monoclonal antibody (1:500).