

## Product datasheet for **TA809772AM**

### PPP3CB Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI2E4]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2E4
Applications:	WB
Recommended Dilution:	WB 1:500~2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 23-96 of human PPP3CB (NP_001135826) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	57.9 kDa
Gene Name:	protein phosphatase 3 catalytic subunit beta
Database Link:	<a href="#">NP_001135826</a> <a href="#">Entrez Gene 19056 Mouse</a> <a href="#">Entrez Gene 24675 Rat</a> <a href="#">Entrez Gene 5532 Human</a> <a href="#">P16298</a>
Synonyms:	CALNA2; CALNB; CNA2; PP2Bbeta
Protein Families:	Druggable Genome, Phosphatase, Transcription Factors




[View online »](#)


**Protein Pathways:**

Alzheimer's disease, Amyotrophic lateral sclerosis (ALS), Apoptosis, Axon guidance, B cell receptor signaling pathway, Calcium signaling pathway, Long-term potentiation, MAPK signaling pathway, Natural killer cell mediated cytotoxicity, Oocyte meiosis, T cell receptor signaling pathway, VEGF signaling pathway, Wnt signaling pathway

**Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY PPP3CB (Cat# [RC227763], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PPP3CB (1:2000) (Cat# [TA809772]). Positive lysates [LY428054] (100ug) and [LC428054] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 8 different cell lines by using anti-PPP3CB monoclonal antibody (1:500).