

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA809927AM

PRMT7 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI10A8]

Product data:

Product Type:	Primary Antibodies	
Clone Name:	OTI10A8	
Applications:	WB	
Recommended Dilution:	WB 1:2000	
Reactivity:	Human, Mouse	
Host:	Mouse	
lsotype:	lgG1	
Clonality:	Monoclonal	
Immunogen:	Full length human recombinant protein of human PRMT7 (NP_061896) produced in HEK293T cell.	
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.	
Concentration:	0.5 mg/ml	
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)	
Conjugation:	Biotin	
Storage:	Store at -20°C as received.	
Stability:	Stable for 12 months from date of receipt.	
Predicted Protein Size:	78.3 kDa	
Gene Name:	protein arginine methyltransferase 7	
Database Link:	<u>NP_061896</u> <u>Entrez Gene 214572 MouseEntrez Gene 54496 Human</u> <u>Q9NVM4</u>	



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	PRMT7 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI10A8] – TA809927AM
--	--

Background:Arginine methylation is an apparently irreversible protein modification catalyzed by arginine
methyltransferases, such as PMT7, using S-adenosylmethionine (AdoMet) as the methyl
donor. Arginine methylation is implicated in signal transduction, RNA transport, and RNA
splicing (Miranda et al., 2004 [PubMed 15044439]). [supplied by OMIM, Mar 2008]. Transcript
Variant: This variant (2) lacks an exon in the 5' UTR and in the 5' coding region compared to
variant 1. The resulting protein (isoform 2) is shorter but has the same N- and C-termini
compared to isoform 1. ##Evidence-Data-START## Transcript exon combination ::
AK304605.1 [ECO:0000332] RNAseq introns :: single sample supports all introns ERS025087,
ERS025093 [ECO:0000348] ##Evidence-Data-END## COMPLETENESS: complete on the 3' end.

Protein Families: Druggable Genome

Product images:

170	-	
130	-	
100	-	
70	-	-
55		
40		
35	-	
25	-	
15	-	
10	-1	

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY PRMT7 (Cat# [RC201672], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PRMT7 (Cat# [TA809927])(1:2000). Positive lysates [LY402728] (100ug) and [LC402728] (20ug) can be purchased separately from OriGene.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US