

Product datasheet for **TA809958AM**


TSHZ1 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI2G9]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2G9
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 29-176 of human TSHZ1(NP_005777) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	112.5 kDa
Gene Name:	teashirt zinc finger homeobox 1
Database Link:	NP_005777 Entrez Gene 110796 Mouse Entrez Gene 307217 Rat Entrez Gene 10194 Human Q6ZSZ6
Synonyms:	CAA; NY-CO-33; SDCCAG33; TSH1
Protein Families:	Transcription Factors



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TSHZ1 ([RC215916], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TSHZ1 (1:2000). Positive lysates [LY417073] (100ug) and [LC417073] (20ug) can be purchased separately from OriGene.