

Product datasheet for TA810066

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ABCD2 Mouse Monoclonal Antibody [Clone ID: OTI3B7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3B7
Applications: WB

Recommended Dilution: WB 1:500

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 418-506 of human

ABCD2 (NP_005155) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 83.1 kDa

Gene Name: ATP binding cassette subfamily D member 2

Database Link: NP 005155

Entrez Gene 26874 MouseEntrez Gene 84356 RatEntrez Gene 225 Human

Q9UBJ2


Synonyms: ABC39; ALDL1; ALDR; ALDRP; hALDR

Protein Families: Druggable Genome
Protein Pathways: ABC transporters





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ABCD2 (Cat# [RC211199], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ABCD2 (Cat# TA810066)(1:500). Positive lysates [LY417474] (100ug) and [LC417474] (20ug) can be purchased separately from OriGene.