

## Product datasheet for **TA810072S**


### **B7H4 (VTCN1) Mouse Monoclonal Antibody [Clone ID: OTI9E9]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI9E9
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 31-117 of human VTCN1 (NP_078902) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	30.7 kDa
Gene Name:	V-set domain containing T cell activation inhibitor 1
Database Link:	<a href="#">NP_078902</a> <a href="#">Entrez Gene 242122 Mouse</a> <a href="#">Entrez Gene 295322 Rat</a> <a href="#">Entrez Gene 79679 Human</a> <a href="#">Q7Z7D3</a>
Synonyms:	B7-H4; B7h.5; B7H4; B7S1; B7X; PRO1291; VCTN1
Protein Families:	Transmembrane



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY VTCN1 (Cat# [RC210360], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-VTCN1 (Cat# [TA810072])(1:2000). Positive lysates [LY403017] (100ug) and [LC403017] (20ug) can be purchased separately from OriGene.