

## Product datasheet for **TA810089S**


### UROD Mouse Monoclonal Antibody [Clone ID: OTI2G7]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2G7
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human UROD (NP_000365) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	40.6 kDa
Gene Name:	uroporphyrinogen decarboxylase
Database Link:	<a href="#">NP_000365</a> <a href="#">Entrez Gene 22275 Mouse</a> <a href="#">Entrez Gene 29421 Rat</a> <a href="#">Entrez Gene 7389 Human</a> <a href="#">P06132</a>
Synonyms:	PCT; UPD
Protein Families:	Druggable Genome
Protein Pathways:	Porphyrin and chlorophyll metabolism



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY UROD (Cat# [RC200529], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UROD (Cat# [TA810089])(1:2000). Positive lysates [LY424758] (100ug) and [LC424758] (20ug) can be purchased separately from OriGene.