

## Product datasheet for **TA810105BM**


### **KRT18P55 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI10E9]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI10E9
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FLJ40504 (NP_775895) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	28.5 kDa
Gene Name:	keratin 18 pseudogene 55
Database Link:	<a href="#">NP_775895</a> <a href="#">Entrez Gene 284085 Human</a>
Synonyms:	MGC138231; MGC138233



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FLJ40504 ([RC211303], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FLJ40504 (1:2000). Positive lysates [LY406513] (100ug) and [LC406513] (20ug) can be purchased separately from OriGene.