

Product datasheet for **TA810106S**


KRT18P55 Mouse Monoclonal Antibody [Clone ID: OTI12H3]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI12H3
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FLJ40504 (NP_775895) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	28.5 kDa
Gene Name:	keratin 18 pseudogene 55
Database Link:	NP_775895 Entrez Gene 284085 Human
Synonyms:	MGC138231; MGC138233



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY FLJ40504 (Cat# [RC211303], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FLJ40504 (Cat# [TA810106])(1:500). Positive lysates [LY406513] (100ug) and [LC406513] (20ug) can be purchased separately from OriGene.