

Product datasheet for **TA810169BM**

Endoglycan (PODXL2) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI7E1]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI7E1
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 366-486 of human PODXL2(NP_056535) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	64.9 kDa
Gene Name:	podocalyxin like 2
Database Link:	NP_056535 Entrez Gene 50512 Human Q9NZ53
Background:	This gene is a member of the CD34 family of cell surface transmembrane proteins, which are characterized by an N-terminal extracellular mucin domain, globular and stalk domains, a single pass transmembrane region, and a charged cytoplasmic tail. The encoded protein is a ligand for vascular selectins. [provided by RefSeq, Oct 2012]
Synonyms:	EG; PODXL2



[View online »](#)

Protein Families: Transmembrane

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY PODXL2 (Cat# [RC208583], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PODXL2 (Cat# [TA810169])(1:2000). Positive lysates [LY414352] (100ug) and [LC414352] (20ug) can be purchased separately from OriGene.