

## Product datasheet for **TA810286BM**


### **DYRK2 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1G8]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI1G8
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 1-210 of human DYRK2 (NP_006473) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	dual specificity tyrosine phosphorylation regulated kinase 2
Database Link:	<a href="#">NP_006473</a> <a href="#">Entrez Gene 8445 Human</a> <a href="#">Q92630</a>
Synonyms:	FLJ21217; FLJ21365
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY DYRK2 (Cat# [RC202864], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DYRK2 (Cat# [TA810286])(1:2000).