

## Product datasheet for **TA810325BM**

### **TROY (TNFRSF19) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1B7]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI1B7
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 30-170 of human TNFRSF19 (NP_061117) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	42.7 kDa
Gene Name:	tumor necrosis factor receptor superfamily member 19
Database Link:	<a href="#">NP_061117</a> <a href="#">Entrez Gene 55504 Human</a> <a href="#">Q9NS68</a>
Synonyms:	TAJ; TAJ-alpha; TRADE; TROY
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction



[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY TNFRSF19 (Cat# [RC214944], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TNFRSF19 (Cat# [TA810325])(1:2000). Positive lysates [LY402704] (100ug) and [LC402704] (20ug) can be purchased separately from OriGene.