

Product datasheet for **TA810435AM**

U2AF1L3 (U2AF1L4) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI8B1]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI8B1
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human U2AF1L4 (NP_659424) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	21.8 kDa
Gene Name:	U2 small nuclear RNA auxiliary factor 1-like 4
Database Link:	NP_659424 Entrez Gene 199746 Human Q8WU68
Synonyms:	U2AF1-RS3; U2AF1L3; U2AF1L3V1; U2AF1RS3; U2af26



[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY U2AF1L4 (Cat# [RC206335], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-U2AF1L4 antibody (Cat# [TA810435]). Positive lysates [LY408126] (100ug) and [LC408126] (20ug) can be purchased separately from OriGene.