

## Product datasheet for **TA810473AM**


### ZNF264 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI10D12]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI10D12
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ZNF264 (NP_003408) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	70.4 kDa
Gene Name:	zinc finger protein 264
Database Link:	<a href="#">NP_003408</a> <a href="#">Entrez Gene 9422 Human</a> <a href="#">O43296</a>
Synonyms:	partial cds; zinc finger protein 264; ZNF264
Protein Families:	Transcription Factors



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZNF264 ([RC214965], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF264 (1:2000). Positive lysates [LY418709] (100ug) and [LC418709] (20ug) can be purchased separately from OriGene.