

## Product datasheet for **TA810476S**


### NCOA7 Mouse Monoclonal Antibody [Clone ID: OTI2E12]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2E12
Applications:	WB
Recommended Dilution:	WB 1:500~2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human NCOA7 (NP_861447) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	106 kDa
Gene Name:	nuclear receptor coactivator 7
Database Link:	<a href="#">NP_861447</a> <a href="#">Entrez Gene 135112 Human</a> <a href="#">Q8NI08</a>
Synonyms:	dj187j11.3; ERAP140; ESNA1; Nbla00052; Nbla10993; NCOA7-AS; TLDC4
Protein Families:	Druggable Genome



[View online »](#)

**Product images:**

Western blot analysis of extracts (35ug) from LOVO cell line by using anti-NCOA7 monoclonal antibody (1:500).