

## Product datasheet for **TA810477**


### **NCOA7 Mouse Monoclonal Antibody [Clone ID: OTI1H1]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	OTI1H1
<b>Applications:</b>	WB
<b>Recommended Dilution:</b>	WB 1:2000
<b>Reactivity:</b>	Human
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG1
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Full length human recombinant protein of human NCOA7 (NP_861447) produced in HEK293T.
<b>Formulation:</b>	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
<b>Concentration:</b>	1 mg/ml
<b>Purification:</b>	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Gene Name:</b>	nuclear receptor coactivator 7
<b>Database Link:</b>	<a href="#">NP_861447</a> <a href="#">Entrez Gene 135112 Human</a> <a href="#">Q8NI08</a>
<b>Synonyms:</b>	dj187j11.3; ERAP140; ESNA1; Nbla00052; Nbla10993; NCOA7-AS; TLDC4
<b>Protein Families:</b>	Druggable Genome



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NCOA7 [RC221673], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NCOA7 (1:2000). Positive lysates [LY405613] (100ug) and [LC405613] (20ug) can be purchased separately from OriGene.