

## Product datasheet for **TA810478M**

### NCOA7 Mouse Monoclonal Antibody [Clone ID: OT1E12]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OT1E12
Applications:	WB
Recommended Dilution:	WB 1:500~2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human NCOA7 (NP_861447) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	nuclear receptor coactivator 7
Database Link:	<a href="#">NP_861447</a> <a href="#">Entrez Gene 135112 Human</a> <a href="#">Q8NI08</a>
Synonyms:	dj187J11.3; ERAP140; ESNA1; Nbla00052; Nbla10993; NCOA7-AS; TLDC4
Protein Families:	Druggable Genome




[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY NCOA7 (Cat# [RC221673], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NCOA7 (Cat# [TA810478])(1:2000). Positive lysates [LY405613] (100ug) and [LC405613] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 3 different cell lines by using anti-NCOA7 monoclonal antibody (1:500).